

Stratify JCV requests are Changing.

Starting in May, Stratify JCV request forms in New Zealand will begin changing, as part of a worldwide update.

The only difference you will see in the collection centres and laboratories is the new style form shown on page 2. It should be marked very clearly at the top “**STRATIFY JCV TEST REQUEST FORM**”

These will still come with 2 barcoded tubes - a primary Red top collection tube and an aliquot tube. Please only use the provided tubes. **Do not label the aliquot tube.**

As they no longer have the collection & processing instructions attached, please see the process as listed below.

Blood Collection

- Collect the barcoded Red top 4mL at the same time as you would any other serum or SST tubes.
- Gently invert the tube 5-10 times.
- When you label the collection tube, do **NOT** cover the 10 digit number or barcode. Do **NOT** label the aliquot tube.
- Please note the time & date of collection on the form, the sampling date listed on the pre-printed form may not be correct.
- Send sample to local lab for processing.

Processing

- Spin primary Red top collection tube for 10 mins, 1800-2000 g.
- Pipette the serum into the krone cryo aliquot tube.
- **Do not label the aliquot tube with Patient information or reference numbers.** The 10 digit number on the tube is adequate identification.
- Store sample in the freezer with its request form.
- Ship the frozen aliquot tube, and the request form together to the Canterbury Health Labs via your normal frozen sample transport method.

Please share this information with all blood collection centres and laboratories collecting & processing samples for Stratify JCV.

If you have any questions around this, please contact us.

Research Team, Canterbury Health Laboratories | Waitaha Canterbury

waea: +64 3 364 1104 pekanga: 81104 | waea pūkoro: +64 027 353 2639
īmēra: chlresearch@cdhb.health.nz

524 Hagley Avenue | PO Box 151, Christchurch 8140



Te Whatu Ora – Health New Zealand

New Zealand

Site No.

YOU SHOULD REGISTER YOUR PATIENT ON-LINE: <https://stratifyjcv.unilabweb.com/>
When using the on-line registration this TRF paperform is not needed - Go paperless Go Green.
All information can be found online including the Sampling handling and courier collection procedure.

1. Patient Identification
Patient identification is mandatory if no full name is provided.

2. Patient last name:
Please write name in full. Use capital letters and the English alphabet

3. Patient first name:
Please write name in full. Use capital letters and the English alphabet

4. Date of Birth:
(e.g. 01-JAN-2021) Enter dat as: DD-MM_YYYY

5. Sex
Male: Female: Other:

6. Which MS medication is the patient currently on?

TYSABRI IV	TYSABRI SC	Aubagio	Avonex	Betaferon	Copaxone	Extavia
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Gilenya	Kesimpta	Lemtrada	Mavenclad	Mayzent	Ocrevus	Plegridy
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Ponvory	Rebif	Tecfidera	Vumerity	Zeposia	Not using MS medication	Other
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

(New to the treatment or treatment break)

7. When was Tysabri initiated?

8. Sampling Date:
E.g. 31-12-1982. Enter date as DD-MM-YYYY.

9. First stratify test?
Yes: No:

Please Note:
The physician must inform the patient that the sample will be sent to Unilabs, Denmark for analysing, and the result is reported to the physician. Upon receipt of the sample, Unilabs will treat patient data confidentially. The sample will be stored by Unilabs for a period of 14 days. After this period, Unilabs reserves the right to destroy the sample without any further notification.

Figure 1: Current Stratify JCV form

Unilabs a.s.
Nygaardsvej 32
2100 Copenhagen
Denmark
Phone: 4533743077
Email: helpdesk@unilabs.com
Stratifyjcv.unilabsweb.com

STRATIFY JCV TEST REQUEST FORM

Carolyn Dutton
Site No.: B1124002
8031

Unilabs Barcode ID: 3773773778 **Reference**

Patient ID: xxx **NHI**

First Name: fake

Last Name: test test

Date of Birth: 14-02-2012

Sex: F

Currently in treatment with: Not currently on medication (new to treatment / treatment break)

Sampling Date: 15-02-2023 **Collection Date**

First Stratify JCV test: Yes

Figure 2: New Stratify JCV Request Form