










OUTTAKES

Tube Guide for Blood Sample Collection

ORDER OF DRAW



| COLOUR | TUBE ADDITIVE | TESTS/DEPT |
|---|-------------------------|----------------------------|
| Blood Cultures | Sterile | Microbiology |
|  | Citrate | INR Coagulation |
|  | CPDA | Tissue typing |
|  | Plain / Clot Activator | |
|  | SST II - Clot Activator | Biochemistry Immunology |
|  | Heparin | Chromosome study |
|  | K2EDTA | Trace metals |
|  | EDTA | Haematology HbA1C |
|  | EDTA | Transfusion Lab |
|  | Fluoride Oxalate | Glucose |